

Improving Flexible Endoscope Reprocessing: the importance of drying and a role for artificial intelligence



Monique Barakat, MD PhD

Associate Professor of Pediatrics & Medicine
Divisions of Pediatric and Adult Gastroenterology and Hepatology
Stanford University School of Medicine

Stanford, CA USA

Intro

This webinar focuses on the issue of endoscope-transmitted infections, highlighting inadequate drying as a contributing factor.¹⁻² We will begin by outlining the challenge of endoscope transmitted infection and identifying key organisms and risk factors, with a particular emphasis on the importance of proper drying.^{1,3-4} Residual fluid is often found in endoscope channels during inspections, underscoring the need for regular evaluations of flexible endoscopes used in clinical settings and more effective approaches to endoscope drying.⁷⁻⁸ We will discuss the

prevalence of inadequate drying and examine practices that contribute to this issue, such as the use of simethicone^{5,15}, endoscope working channel⁷⁻⁸ damage and suboptimal eprocessing methods.^{3,9}

We will compare various drying techniques - including drying cabinets¹⁻², automated systems, and manual methods - detailing their respective advantages and disadvantages while also presenting strategies to minimize human error and enhance patient safety.^{3,13-14}


Finally, the session will address the practical challenges of endoscope inspection and introduce the emerging role of AI-assisted inspection, supported by pilot data and demonstration videos.¹⁶



ASP Advanced Sterilization Products

Endoscope transmitted infections

Endoscope-transmitted infections (ETIs) pose a significant threat to patient safety, as they can lead to severe health complications, including sepsis and prolonged hospital stays.² The complexity of endoscope design makes thorough cleaning and disinfection challenging, increasing the risk of residual pathogens and endoscope transmitted infections.^{1,3} Inadequate adherence to reprocessing protocols can result in the transmission of harmful bacteria, viruses, and fungi, which may originate from patient flora or environmental sources.^{2,5} Outbreaks of ETIs have raised awareness and highlighted the need for stringent infection control measures in healthcare settings.^{2,13} Addressing the risks associated with ETIs is crucial for improving patient outcomes and ensuring the safety of endoscopic procedures.^{3,13}



Endoscope Reprocessing is a complex process

- 1 Leak tests
- 2 Manual cleaning
- 3 Automated cleaning
- 4 High Level Disinfection / Liquid Chemical Sterilisation
- 5 Transport and or Storage & Drying
- 6 Clinical Application/ use

“The complexity of endoscope design makes thorough cleaning and disinfection challenging, increasing the risk of residual pathogens and endoscope transmitted infections.^{1,3}”

Addressing the risks associated with ETIs is crucial for improving patient outcomes and ensuring the safety of endoscopic procedures.^{3,13}



The importance of drying endoscopes

Complete drying of endoscopes is essential for preventing endoscope-transmitted infections (ETIs), because residual moisture can promote microbial growth and biofilm formation within endoscopes.^{2,4,10,12}

Inadequate drying increases the likelihood of pathogens surviving the disinfection process, thereby heightening the risk of infection transmission to patients.^{2,3,5,13}

Prioritizing endoscope drying during reprocessing and implementing effective drying protocols is essential to enhance the overall safety of endoscopic procedures and reduce the potential for healthcare-associated infections.^{1,3,10,12}

Prioritizing endoscope drying during reprocessing and implementing effective drying protocols is essential to enhance the overall safety of endoscopic procedures and reduce the potential for healthcare-associated infections.^{1,3,10,12}




“Complete drying of endoscopes is essential for preventing endoscope-transmitted infections (ETIs), because residual moisture can promote microbial growth and biofilm formation within endoscopes.^{2,4,10,12}”



Image source: St George Endoscopy Department - drying cabinet - courtesy Mr David Bellamy, Sydney, Australia, 2025 Google photo search "flexible endoscope drying cabinet" access on 20250522 <https://www.aorn.org/outpatient-surgery/article/2022-September-endoscope-cabinets> Source: ISO 15883-4:2008 - Washer-disinfectors for thermolabile endoscopes <https://www.iso.org/standard/63696.html> / GENCA Guidelines 2025 (Australia), https://www.genca.org/public/5/files/Nurses%20info/IPCE%202021_Feb2022update.pdf / AAMI ST91:2021 (USA), <https://www.aami.org/ST91> / AS 5369:2023 (Australia) - Reprocessing of reusable medical devices in health service organizations <https://www.standards.org.au/blog/spotlight-on-as-5369-2023>

Endoscope channel inspection has revealed that drying is often inadequate

Inspections of endoscope channels have consistently revealed that a significant proportion of endoscopes are inadequately dried.^{7,8,16} This inadequate drying poses a serious risk for endoscope-transmitted infections.^{2,3,5,13} The presence of residual moisture in these channels can facilitate the growth of harmful pathogens, undermining the effectiveness of high-level disinfection protocols.^{4,10,12,15} Endoscope working channel damage may predispose to fluid retention, with areas of damage serving as a sanctuary for fluid and pathogens.^{7,8,16} Recognizing and addressing the prevalence of inadequate drying is crucial for enhancing patient safety and preventing healthcare-associated infections during endoscopic procedures.^{1,3,5,10}



Source: Barakat MT, Girotra M, Huang RJ, Banerjee S. Scoping the scope: endoscopic evaluation of endoscope working channels with a new high-resolution inspection endoscope. GIE 2018 Oct;88(4):601-611.e1. PMID: 29425885

Spectrum of Endoscope Damage

Superficial scratches

Adherent Peel

Deep scratches

Burns

Channel Buckling


Stains

Perforation

“Inspections of endoscope channels have consistently revealed that a significant proportion of endoscopes are inadequately dried.”^{7,8,16}


Leak Test Failure

Scope Inspection Results: Endoscope Damage



- Superficial scratches (98.5%)
- Scratches with adherent peel (76.5%).
- Rare channel buckling
- No channel perforations, stains, or burns.

Scope Inspection Results: Fluid & Debris



- Minor punctate debris was common (>90%)
- Fluid droplets were noted in 42.6% of endoscopes after reprocessing and drying.

Recognizing and addressing the prevalence of inadequate drying is crucial for enhancing patient safety and preventing healthcare-associated infections during endoscopic procedures.^{1,3,5,10}



Topic • 04

Evidence-based approaches to effective endoscope drying

Identification of optimal approaches for endoscope drying has been an area of active research and evidence-based approaches for effective drying of endoscopes during reprocessing have been identified.^{1,3,10,12} These approaches include implementing fixed duration of manual or automated forced air drying systems, which can significantly reduce drying time and enhance the efficacy of the drying process by delivering high-velocity air into the endoscope channels.^{10,12}

Qualitative Visual Assessment of Residual Fluid after Automated Drying



Manual Forced Air Drying



Automated Drying 5 min



Automated Drying 10 min



Source: Barakat M T, Huang R J, Banerjee S. Comparison of automated and manual drying in the elimination of residual endoscope working channel fluid after reprocessing (with video) *Gastrointest Endosc.* 2019;89:124–132.



“Identification of optimal approaches for endoscope drying has been an area of active research and evidence-based approaches for effective drying of endoscopes during reprocessing have been identified.^{1,3,10,12}”

The use of automated drying cabinets that provide controlled environments to ensure thorough moisture removal has also been validated.¹² With any protocol implemented, regular monitoring and validation of drying protocols is essential to maintain compliance and ensure optimal drying outcomes, ultimately minimizing the risk of endoscope-transmitted infections.^{1,3,11}

Automated Drying Cabinets Effectively Dry Endoscopes & Bronchoscopes

Fluid from
Scopes, if
present, on
Cobalt
Chloride Paper



No Fluid



Fluid Present



Cobalt
Chloride Paper
Analysis



Dry



Wet

		Cabinet													
				Standard	Automated	Standard	Automated	Standard	Automated	Standard	Automated	Standard	Automated	Standard	Automated
Bronchoscope (BF-3C20)	Internal Channels														
	External Channels														
Duodenoscope (TJF-160F)	Internal Channels														
	External Channels														
Colonoscope (CF-Q160AL)	Internal Channels														
	External Channels														
				0.5h		1h		2h		3h		24h			

Source: Perumpail RB, Marya NB, McGinty BL, Muthusamy VR. Endoscope reprocessing: Comparison of drying effectiveness and microbial levels with an automated drying and storage cabinet with forced filtered air and a standard storage cabinet. Am J Infect Control. 2019 Sep;47(9):1083-1089. doi: 10.1016/j.ajic.2019.02.016. Epub 2019 Apr 6. PMID: 30962022.

With any protocol implemented, regular monitoring and validation of drying protocols is essential to maintain compliance and ensure optimal drying outcomes, ultimately minimizing the risk of endoscope-transmitted infections.^{1,3,11}




Topic • 05

Role for AI in endoscope reprocessing

In conclusion, addressing endoscope-transmitted infections (ETIs) is critical for ensuring patient safety, as these infections can lead to severe health complications and highlight the need for stringent infection control measures.^{2,3,5,13} Complete drying of endoscopes is essential to prevent microbial growth and biofilm formation, as inadequate drying significantly increases the risk of infection transmission.^{2,4,10,12,15} Inspections have revealed that many endoscopes are inadequately dried, underscoring the importance of implementing effective drying protocols and regular monitoring of moisture elimination within the endoscope channels to enhance patient safety.^{7,8,10,12,16} Finally, the integration of artificial intelligence (AI) in the inspection process offers a promising avenue for improving efficiency and accuracy, ultimately contributing to better outcomes in endoscope reprocessing and the safety of patients who undergo endoscopy.¹⁶



The inspection of endoscope working channels in a clinical setting presents significant demands in terms of time to conduct the inspections and expertise to identify, quantify and summarize working channel findings.^{7,8,16}



Source: Barakat MT, Girotra M, Huang RJ, Banerjee S. Scoping the scope: endoscopic evaluation of endoscope working channels with a new high-resolution inspection endoscope. GIE 2018 Oct;88(4):601-611.e1. PMID: 29425885



“To streamline this process of endoscope working channel inspections and enhance accuracy, the integration of artificial intelligence (AI) is emerging as a promising solution, offering the potential to automate and improve inspection efficiency.”¹⁶

Conclusion

In conclusion, addressing endoscope-transmitted infections (ETIs) is critical for ensuring patient safety, as these infections can lead to severe health complications and highlight the need for stringent infection control measures.^{2,3,5,13} Complete drying of endoscopes is essential to prevent microbial growth and biofilm formation, as inadequate drying significantly increases the risk of infection transmission.^{2,4,10,12,15} Inspections have revealed that many endoscopes are inadequately dried, underscoring the importance of implementing effective drying protocols and regular monitoring of moisture elimination within the endoscope channels to enhance patient safety.^{7,8,10,12,16} Finally, the integration of artificial intelligence (AI) in the inspection process offers a promising avenue for improving efficiency and accuracy, ultimately contributing to better outcomes in endoscope reprocessing and the safety of patients who undergo endoscopy.¹⁶

Take home messages



1. Endoscope-transmitted infections (ETIs) pose a significant threat to patient safety, necessitating stringent infection control measures to prevent severe health complications.^{2,3,5,13}

2. Complete drying of endoscopes is essential to prevent microbial growth and biofilm formation, thereby reducing the risk of endoscope-transmitted infections (ETIs).^{2,4,10,12,15}

3. Inspections have consistently shown that inadequate drying of endoscopes significantly increases the risk of endoscope-transmitted infections, highlighting the need for improved drying protocols.^{7,8,16}

4. Implementing evidence-based drying approaches, such as automated drying systems and regular monitoring, is crucial for ensuring optimal drying outcomes and minimizing the risk of ETIs.^{1,3,10,11,12}

5. The integration of artificial intelligence (AI) in endoscope inspections offers a promising solution to enhance efficiency and accuracy, addressing the significant demands of the inspection process.¹⁶

20
25

ASP SUMMIT
CONTINUOUS EDUCATION



Sources:

1. Barakat MT, Banerjee S. Novel Algorithms for Reprocessing, Drying and Storing Endoscopes. *Gastrointest Endosc Clin N Am*. 2020 Oct;30(4):677-691. PMID: 32891225.
2. Kovaleva J, Peters FT, van der Mei HC et al. Transmission of infection by flexible gastrointestinal endoscopy and bronchoscopy. *Clin Microbiol Rev*. 2013;26:231-254.
3. Day LW, Muthusamy VR, Collins J, Kushnir VM, Sawhney MS, Thosani NC, Wani S. Multisociety guideline on reprocessing flexible GI endoscopes and accessories. *Gastrointest Endosc*. 2021 Jan;93(1):11-33.e6. doi: 10.1016/j.gie.2020.09.048. Erratum in: *Gastrointest Endosc*. 2021 Nov;94(5):1019. doi: 10.1016/j.gie.2021.09.026. PMID: 33353611.
4. Alfa MJ, Sitter DL. In-hospital evaluation of contamination of duodenoscopes: a quantitative assessment of the effect of drying. *J Hosp Infect* 1991; 19: pp. 89-98.
5. Cheung DY et al. Multidisciplinary and Multisociety Practice Guideline on Reprocessing Flexible Gastrointestinal Endoscopes and Endoscopic Accessories. *Clin Endosc*. 2020 May;53(3):276-285. doi: 10.5946/ce.2020.106. Epub 2020 May 29. PMID: 32506894.
6. Higa JT, Choe J, Toms D et al. Optimizing duodenoscope reprocessing: rigorous assessment of a culture and quarantine protocol. *Gastrointest Endosc*. 2018;88:223-229.
7. Barakat MT, Girotra M, Huang RJ, Banerjee S. Scoping the scope: endoscopic evaluation of endoscope working channels with a new high-resolution inspection endoscope. *GIE* 2018 Oct;88(4):601-611.e1. PMID: 29425885
8. Thaker A M, Kim S et al. Inspection of endoscope instrument channels after reprocessing using a prototype borescope. *Gastrointest Endosc*. 2018;88:612-619.
9. Day L W, Muthusamy V R, Collins J et al. Multisociety guideline on reprocessing flexible GI endoscopes and accessories. *Gastrointest Endosc*. 2021;93:11-33 e16.
10. Barakat M T, Huang R J, Banerjee S. Comparison of automated and manual drying in the elimination of residual endoscope working channel fluid after reprocessing (with video) *Gastrointest Endosc*. 2019;89:124-132.
11. ASGE Technology Committee; Komanduri S, Abu Dayyeh BK, Bhat YM, Chauhan SS, Gottlieb KT, Hwang JH, Konda V, Lo SK, Manfredi M, Maple JT, Murad FM, Siddiqui UD, Wallace MB, Banerjee S. Technologies for monitoring the quality of endoscope reprocessing. *Gastrointest Endosc*. 2014 Sep;80(3):369-73. doi: 10.1016/j.gie.2014.01.044. Epub 2014 Jul 8. PMID: 25012558.
12. Perumpail RB, Marya NB, McGinty BL, Muthusamy VR. Endoscope reprocessing: Comparison of drying effectiveness and microbial levels with an automated drying and storage cabinet with forced filtered air and a standard storage cabinet. *Am J Infect Control*. 2019 Sep;47(9):1083-1089. doi: 10.1016/j.ajic.2019.02.016. Epub 2019 Apr 6. PMID: 30962022
13. Beilenhoff U, Biering H, Blum R et al. Reprocessing of flexible endoscopes and endoscopic accessories used in gastrointestinal endoscopy: Position Statement of the European Society of Gastrointestinal Endoscopy (ESGE) and European Society of Gastroenterology Nurses and Associates (ESGENA) – Update 2018. *Endoscopy*. 2018;50:1205-1234.
14. Day LW, Muthusamy VR, Collins J, Kushnir VM, Sawhney MS, Thosani NC, Wani S. Multisociety guideline on reprocessing flexible GI endoscopes and accessories. *Gastrointest Endosc*. 2021 Jan;93(1):11-33.e6. doi: 10.1016/j.gie.2020.09.048. Erratum in: *Gastrointest Endosc*. 2021 Nov;94(5):1019. doi: 10.1016/j.gie.2021.09.026. PMID: 33353611.
15. Barakat M T, Huang R J, Banerjee S. Simethicone is retained in endoscopes despite reprocessing: impact of its use on working channel fluid retention and adenosine triphosphate bioluminescence values (with video) *Gastrointest Endosc*. 2019;89:115-123
16. Barakat MT, Girotra M, Banerjee S. Initial application of deep learning to borescope detection of endoscope working channel damage and residue. *Endosc Int Open*. 2022 Jan 14;10(1):E112-E118. doi: 10.1055/a-1591-0258. PMID: 35047341

ASP Advanced
Sterilization
Products

ASP International GmbH, Im Majorenacker 10,
Schaffhausen Switzerland
©ASP 2025. All Rights Reserved.

asp.com



SCAN ME

Advanced Sterilization Products (ASP) is not responsible nor can be held liable for the accuracy of the information and data provided by the health care professionals (HCPs) who have presented and/or published such information and data that ASP and permitted thereafter to ASP to incorporate into this work product. It remains the HCPs responsibility to ensure that their presentation and/or publications are supported by factual and sourced evidences, in light of the best scientific knowledge and experience available at the time of the presentation and/or publication. Capitalized product names are trademarks of ASP Global Manufacturing GmbH.

SM-2500031-01-1